

MATLAB code documentation: “Penumbra quantification from MR SWI-DWI mismatch and its comparison with MR ASL PWI-DWI mismatch in patients with acute ischemic stroke”



IIT Delhi

Indian Institute of Technology Delhi

Code prepared by:

Rupsa Bhattacharjee, Dr. Anup Singh, et.al.

MedImg Group,

Center for Biomedical Engineering (CBME) IIT Delhi

This code sharing module is prepared for Research usage strictly




FIGURE 1 Flowchart of proposed methodology for penumbra quantification from MR SWI-DWI mismatch and its comparison with MR ASL-PWI-DWI mismatch

Block A: Preprocessing: Read the data, resize, register (Line: 3-255) : Descriptions and examples given for small segments of the block

1st Segment (Line 6 -24)

Description: (Alternatively, use lab-specific data reading practices)

- Insert the folder location where the DICOM data is stored.
- Enter the total number of slices/images for all the four sequences: SWI, pCASL (PWI/CBF), ADC and FLAIR.
- Select all the DICOM Images from each of the respective folders.
- Define the resolution of SWI. We have worked with 1008 X 1008.


```
%% Informations about the Input Data

InputDir = uigetdir('Select patient data folder'); %% Insert the folder location, where the DICOM data is kept
cd(InputDir)

SliceTotalSWI = input('Total Number of Slices in SWI = '); %%Input the total number of slices/images in SWI sequence, example n = 130/140/150
SliceTotalpCASL = input('Total Number of Slices in pCASL = '); %%Input the total number of slices/images in pCASL sequence, example n = 14
SliceTotalDWI = input('Total Number of Slices in ADC = '); %%Input the total number of slices/images in ADC sequence, example n = 25
SliceTotalFLAIR = input('Total Number of Slices in Flair = '); %%Input the total number of slices/images in FLAIR sequence, example n = 25


%%Select all the dicom images from each of the folders
[filenameSWI_mag, pathnameSWI_mag] = uigetfile('multiselect', 'on', '*.dcm', 'select SWI mag images');
[filenameCASL, pathnameCASL] = uigetfile('multiselect', 'on', '*.dcm', 'select 3D ASL images');
[filenameADC, pathnameADC] = uigetfile('multiselect', 'on', '*.dcm', 'select ADC images');
[filenameFLAIR, pathnameFLAIR] = uigetfile('multiselect', 'on', '*.dcm', 'select FLAIR images');

%%Resolutions of SWI images (In our case this is 1008 X 1008), can be
%%modified as per input data.
H = 1008; H_SWI = 1008; W = 1008; W_SWI = 1008;
```



Total Number of Slices in SWI = 145
 Total Number of Slices in pCASL = 14
 Total Number of Slices in ADC = 25
 Total Number of Slices in Flair = 25

Enter in Command window



Use Ctrl + A and select all the slices/images in each of the four sequence folders

2nd Segment (Line 25 – 83)

Description:

- i) Read the input data from location previously specified.
- ii) Run either of the modules depending on whether the data is in Enhanced DICOM Format or Classic DICOM format
- iii) Define the Starting and ending slice range: where the diffusion restriction is seen on ADC (Ex: 10 to 19 in given example)
- iv) Define the Starting and ending slice range: where the Hypo-perfusion is seen on ASL_CBF
- v) Define the Starting and ending slice range: where the PHVS is seen on SWI
- vi) Define the Starting and ending slice range: where the Abnormality is seen on FLAIR

```

%% READ the Input Data
%%IF Dicom datasets are ENHANCED DICOM: Run this module
SI_SWI_mag1 = double(dicomread(fullfile(pathnameSWI_mag,filenameSWI_mag)));
SI_pCASL1 = double(dicomread(fullfile(pathnamepCASL,filenamepCASL)));
SI_ADC1 = double(dicomread(fullfile(pathnameADC,filenameADC)));
SI_FLAIR1 = double(dicomread(fullfile(pathnameFLAIR,filenameFLAIR)));
for i=1:(SliceTotalSWI)
    SI_SWI_mag(:,i) = SI_SWI_mag1(:,i);
end
for i=1:SliceTotalpCASL
    SI_pCASL(:,i) = SI_pCASL1(:,i);
end
for i=1:SliceTotalDWI
    SI_ADC(:,i) = SI_ADC1(:,i);
end
for i=1:SliceTotalFLAIR
    SI_FLAIR(:,i) = SI_FLAIR1(:,i);
end

```


For Enhanced DICOM

```

%%IF Dicom datasets are CLASSIC DICOM: Run this module
for i=1:SliceTotalSWI
    SI_SWI_mag(:,i) = double(dicomread(fullfile(pathnameSWI_mag,filenameSWI_mag(i))));
end
for i=1:SliceTotalpCASL
    SI_pCASL(:,i) = double(dicomread(fullfile(pathnamepCASL,filenamepCASL(i))));
end
for i=1:SliceTotalDWI
    SI_ADC(:,i) = double(dicomread(fullfile(pathnameADC,filenameADC(i))));
end
for i=1:SliceTotalFLAIR
    SI_FLAIR(:,i) = double(dicomread(fullfile(pathnameFLAIR,filenameFLAIR(i))));
end

```

For Classic DICOM



Example : Slice number 10 to 19 shows Diffusion restriction

3rd Segment (Line 84 – 255)

Description: Run the modules/lines, no inputs required, change values as necessary

- i) Resize ADC , FLAIR & 3D pCASL images as per SWI resolution (Highest matrix size)
- ii) Create sub-matrix of ADC FLAIR & 3D pCASL images, containing only relevant abnormal slices (Rest of them are not taken into consideration)
- iii) Create MIP (minimum intensity projection for SWI) for every 6 slices. In this code, SWI is of 1 mm slice thickness, rest sequences (ADC, FLAIR, pCASL) are of 6 mm slice thickness. Hence

the MIP is created for every 6 SWI slices, to match the rest. This can be modified as per user-specific sequence parameters.


- iv) Registration of ADC, FLAIR and ASL with SWI MIP

Block B: Enhancement of PHVS Features (Line: 256-306)

Description: Run the modules/lines, no inputs required, change values if necessary

- i) Frangi 2D Filter implementation on SWI MIP
- ii) Frangi 3D implementation on SWI MIP
- iii) Frangi 2D and 3D Combination and Top Hat Morphological transformation
- iv) Gaussian Filter + Morphological Dilation
- v) Final PHVS feature extracted image is generated

PHVS feature extracted





Block C: Penumbra Probability mask creation and automatic segmentation of SWI Penumbra (Line: 309-368)

Description: Run the modules/lines, change values if necessary

- i) Mark the Stroke occurring hemisphere (SOH) & Contralateral Normal appearing hemisphere (CNAH) on SWI MIP images (Each MIP Slice)
- ii) Multiplying the ROI's to PHVS feature images

- iii) ADC Core segmentation: We have automatized upon human threshold, observed, and reported. Can be manually segmented or changed as per animal models. This core mask works as a seed point.
- iv) Binarize, vessel density maps and automatize using functions
- v) **SWI_penumbra is the final segmented BW mask.**



SWI_Penumbra mask, overlaid on SWI MIP

Block D: Volumes and Mismatch ratios calculation (Line: 369-461)

Description: Run the modules/lines, change values if necessary

- Final output values are **Volumes** : ICV_in_ml, PPV_in_ml, SPV_in_ml, FLV_in_ml, **Ratios**: Mismatch_Ratio_SWI_ADC, Mismatch_Ratio_pCASL_ADC, Mismatch_Ratio_FLAIR_ADC

- For our study, FOV was fixed at 230 X 230. If the FOVs are changed for any sequence, accordingly those values need to be updated in this section.
- For our study, slice thickness and gaps were fixed at 6 mm and 0 mm respectively. If these values are changed for any sequence, accordingly those values need to be updated in this section.

If you use the framework, please cite the following Reference:

[1] R. Bhattacharjee, R. K. Gupta, B. Das, V. K. Dixit, P. Gupta, and A. Singh, “Penumbra quantification from MR SWI-DWI mismatch and its comparison with MR ASL PWI-DWI mismatch in patients with acute ischemic stroke,” *NMR Biomed.*, no. December 2020, pp. 1–15, 2021, doi: 10.1002/nbm.4526.

Contact Details: For any further clarifications, please contact:

- Rupsa Bhattacharjee, rupsa.bhattacharjee1@gmail.com, +91-8800172999
- Dr. Anup Singh, Anup.Singh@cbme.iitd.ac.in, +91-7042311855